

	N Class non-vacuum		S Class vacuum		B Class vacuum		
	Solid loads Unwrapped		Solid and Hollow Loads S/M Weight Porous Loads Wrapped or unwrapped		Solid and Hollow Loads S/M/L Weight Porous Loads Wrapped or unwrapped		
Chamber Size	18 Ltr	23 Ltr	18 Ltr	23 Ltr	18 Ltr	23 Ltr	60 Ltr
Chamber Length	38 cm	53 cm	38 cm	53 cm	38 cm	53 cm	60 cm
Tray Dimensions Total number of trays Size (mm) Size (mm)	4 or 5 185x285	4 or 5 185x440	4 or 5 185x285	4 or 5 185x440	4 or 5 185x285	4 or 5 185x440	4 2 @ 310x550 2 @ 220x550
Sterilizable Weight Unwrapped (kg) Wrapped (kg) Porous (kg)	3 - -	4.5 - -	2.5 1.5 0.5	3.5 2 0.5	4 2.5 1	4.5 3.5 1.5	12 12 6
Outer Dimensions Depth (mm) Width (mm) Height (mm) Weight (kg)	615 505 400 43	690 505 400 48	615 505 400 45	690 505 400 50	615 505 400 45	690 505 400 50	783 680 1600 205